

Minidop ES-100VX

POCKET DOPPLER WITH GREAT SOUNDS



- **Pocketable, Hand-Held Size**
- **Wide probe selection for user flexibility**
Clinics: 2, 4, 5, 8 and 10MHz
Surgical*: 8, 10 and 20MHz

- **Great sensitivity**
- **Convenient probe activation button**
for turning the unit ON and OFF

CLINICAL APPLICATIONS

- Blood flow detection

ACCESSORIES

- Carrying case
- Battery
- Ultrasonic gel

DOPPLER PROBES

Standard (Clinic): 2, 4, 5, 8 and 10MHz



Optional probes:

- Pencil: 8MHz
- Flat: 8MHz

Sterilizable probes*: (Amplifier required)

- Pencil: 8, 10 and 20MHz
- Autoclavable: 8MHz

Other probes available

**All Sterilizable probes are not available in European Union Countries*

SPECIFICATIONS

Power Supply: DC 9V, alkaline battery
Battery life: Approx. 3 hours
Automatic shut off: No signal: 3 min.
 No signal: 3min.
Speaker output: Max. 300mW or more
Multi frequency: 2, 4, 5, 8, 10 and 20MHz
External outputs: Headset connector

Dimensions: 78(W) x 141(L) x 27(H) mm
 (Not including probe holder)
Weight: Approx. 240g
 (Including battery and probe)
Storage and Transport environment:
 0 to 50°C 85% humidity or less with
 no condensation
Operating environment: 10 to 37°C
Electrical safety: Conform to IEC60601-1

Specifications are subject to change.



Hadeco, Inc.
 2-7-11, Arima, Miyamae-ku, Kawasaki,
 Kanagawa, 216-0003, Japan
 TEL: +81 44 877 4361
 FAX: +81 44 855 7301
<http://www.hadeco.co.jp>



080-00125-09
 Printed in Japan